

Dip & Read Tests

For Schools

Rapid Tests for Schools

- Economical
- High educational value
- Easy to use



MACHERY-NAGEL

www.mn-net.com







Since 1911

MN Dip & Read Tests for Schools

Test strips motivate kids to discover their environment in many different ways. They can easily check the water quality in the local river, the ozone concentration in the summer air, the vitamin C content of fruits etc. The easy Dip & Read procedure of test strips safely provides reliable results in a few seconds, making MN products the perfect choice for an exiting but safe educational environment.



Product information

Test	Product	Application	Gradation	REF
	Ozone Test strip	<ul style="list-style-type: none"> Ozone is formed on sunny days from oxygen and nitrogen oxides Easy test to measure the local ozone concentration Provides results in 10 minutes 	< 90 · 90 – 150 · 150 – 210 · > 210 µg/m ³ O ₃	907 36
	AQUADUR® 3 - 25	<ul style="list-style-type: none"> Easy determination of water hardness 5 different pads for easy read off Provides valuable information about the local water 	<3 · >5 · >10 · >15 · >20 · >25 °d	912 01
	QUANTOFIX® Nitrate/Nitrite	<ul style="list-style-type: none"> Fast and easy detection of nitrate and nitrite Check water quality in the local river Detects nitrate in vegetables and potatoes 	0 · 10 · 25 · 50 · 100 · 250 · 500 mg/l NO ₃ ⁻ 0 · 1 · 5 · 10 · 20 · 40 · 80 mg/l NO ₂ ⁻	913 13
	pH-Fix 0-14	<ul style="list-style-type: none"> Easy detection of pH value in surface water, rivers, fruit juice etc. Brilliant colour charts for clear results Plastic strip prevents skin contact with the sample 	pH 0 · 1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	921 10
Others				
	Litmus paper, neutral	Basic determination of acids and bases		911 07
	Potassium iodide starch paper	Detection of oxidisers like free chlorine and nitrite		907 54
	QUANTOFIX® Ascorbic acid	For vitamin C in juice or fruits		913 14
	Oil test paper	Oil detection in water		907 60
	Sulfide test paper	Detection of hydrogen sulfide		907 61
	Medi-Test Glucose	Specific test for glucose		930 24
	Watesmo	Easy detection of water		906 09

Contact

Your local distributor




Image Credits © Robt. BlueOrange Studio, Noam - Fotolia.com
KATEN100099/Schools en1/2/0/10.2009 PD
Printed in Germany

MACHEREY-NAGEL

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany
 Germany and international: Tel.: +49 (0) 24 21 96 90 Fax: +49 (0) 24 21 96 91 99 e-mail: sales-de@mn-net.com
 Switzerland: Tel.: +41 (0) 62 388 55 00 Fax: +41 (0) 62 388 55 05 e-mail: sales-ch@mn-net.com
 France: Tel.: +33 (0) 3 88 68 22 68 Fax: +33 (0) 3 88 51 76 88 e-mail: sales-fr@mn-net.com
 USA: Tel.: +1 484 821 0984 Fax: +1 484 821 1272 e-mail: sales-us@mn-net.com

